



Stereotests

Stereo depth perception tests are an effective and easy to use method of screening stereoscopic vision for all ages. They help to detect visual problems such as stereopsis, amblyopia, suppression and strabismus.



THE ORIGINAL STEREO FLY

(GOLD STANDARD FOR STEREO DEPTH PERCEPTION TESTING)

- Evaluate both gross stereopsis (355 to 700 seconds of arc) and fine depth perception
- Graded circle test (800 to 40 seconds of arc)
- Animal test for children (400 to 100 seconds of arc)

Test includes:

- 1 pair of Standard 3-D Viewers
- Answer key on back cover




STEREO BUTTERFLY

- Tests and evaluates random dot stereo depth perception for both gross stereopsis and fine depth perception
- Graded circle test (800 to 40 seconds of arc)
- Animal test for children (400 to 100 seconds of arc)

Test includes:

- 1 pair of Standard 3-D Viewers
- Answer key on back cover



RANDOT®

- Tests and evaluates stereo depth perception for both gross stereopsis evaluation (500 & 250 seconds of arc) using six geometric forms from random dot backgrounds
- Graded circle test (400 to 20 seconds of arc)

Test includes:

- 1 pair of Standard 3-D Viewers
- Answer key on back cover