



## SPECIFICATION OF TRIAL BOX WITH TRIAL FRAME

### TRIAL BOX (WOODEN)

Trial box contains trial lenses arranged in pairs with prisms and accessories.

### TRIAL LENSES

Reduced aperture lenses of 20 mm diameter mounted in the metal rim of 38mm diameter, Plano convex and Plano concave.

- 1) Spherical lenses : 32 pairs present in plus and minus power. Ranges :  $\pm 0.12D$  to  $\pm 20.00D$ .  
 Variations: 0.25 to 4.00 D in increment of 0.25 D.  
               4.00D to 6.00 D in increment of 0.50 D.  
               6.00 to 14.00 D in increment of 1.00 D.  
               14.00 to 20.00 D in increment of 2.00D.
- 2) Cylindrical lenses : 19 pairs present in plus and minus power. Ranges :  $\pm 0.25D$  to  $\pm 6.00D$ .  
 Variations: 0.25 to 3.50D in increment of 0.25D.  
               3.50 to 6.00 in increment of 0.50D.
- 3) Prisms power range of 0.5 D, 1 to 6.00D in 1D steps, 6 to 12.00D in 2D steps.

### ACCESSORIES:

1. Trial frame